


2 mL Stainless Steel Loading Guidelines

Bead Ruptor Elite™ - Maximum Speed: 6 m/s
Maximum of 8 tubes can be loaded.


Tubes must be balanced in the tube carriage.



Place tubes spaced evenly
In the tube carriage.

Bead Ruptor™ 96 - Maximum Speed: 30 Hz
Minimum of 4 tubes must be loaded.

Tubes must be balanced in the Tube Holder.



Place tubes spaced evenly
In the Tube Holder.



2 mL Stainless Steel Loading Guidelines

WARNING:

- Not Compatible with Bead Ruptor™ 4 Homogenizer or Bead Ruptor™ 12 Homogenizer
- Cleaning – Autoclavable, wash out with water / 70% Ethanol until clean. Avoid long contact time with bleach
- Metal Lysing tube containing a 6 x 12.5 mm Stainless Steel Slug
- Changing the tube fill from the included 6 x 12.5 mm Stainless Steel Slug to any other fill is at the user's own risk and responsibility.

