

# Nucleic Acid Extraction from HepG2 Cells: An Integrated Solution Utilizing the Bead Ruptor Elite™ Homogenizer, chemagic™ 360 Nucleic Acid Extractor and JANUS® G3 Automated Workstation

Gabriella Ryan, M.S.; Rodney Nash, Ph.D  
PerkinElmer

Bead Ruptor Elite™ Bead Mill Homogenizer



## Summary

Applications featuring Next Generation Sequencing, quantitative PCR, or other genomics-based research using *in-vitro* cultured cells requires a nucleic acid extraction upstream of aforementioned methods. The Bead Ruptor Elite™ Bead Mill Homogenizer and 2 mL bead kits are a streamlined sample preparation solution that provides efficiency in cell lysis, releasing nucleic acid for extraction on the chemagic™ 360 automated nucleic acid extractor.

Herein, we outline integration of the Bead Ruptor Elite™ Bead Mill Homogenizer with the chemagic™ 360 nucleic acid extractor and JANUS® G3 automated workstation, creating an end-to-end solution for sample preparation of *in-vitro* cultured cells and nucleic acid extraction workflows.

## Materials and Methods

### Equipment

- Bead Ruptor Elite™ Bead Mill Homogenizer (Cat # 19-042E)
- Bead Ruptor Elite™ Bead Mill Homogenizer 2 mL Tube Carriage (Cat # 19-373)
- Bead Ruptor Elite™ Bead Mill Homogenizer Finger Plate (Cat # 19-370)
- Bead Ruptor Elite™ Bead Mill Homogenizer 48 position 2 mL Tube Carriage (Cat # 19-378)
- Bead Ruptor Elite™ Bead Mill Homogenizer 48 position 2 mL Tube Carriage Finger Plate (Cat # 19-370-248)
- Hard Tissue Homogenizing Mix 2.8 mm Ceramic Beads (Cat # 19-628)

# Materials and Methods

## HepG2 Cell Culture and Cell Counting


HepG2 cells (ATCC, Cat # HB-8065) were grown in DMEM (Gibco, Cat # 11965118) supplemented with 7% FBS (Gemini Bio, Cat # 900-108), 1% L-glutamine (Gemini Bio, Cat # 400-106) and 1% Penicillin:Streptomycin (Gemini Bio, Cat # 400-109) on tissue culture treated 75 cm<sup>2</sup> flasks (Corning, Cat # 430641U). Cells were grown to confluency at 37°C supplemented with 5% CO<sub>2</sub>, and were detached from the plate using Accutase (Gemini Bio, Cat # 400-158). HepG2 cells were counted using the Cell Size assay on the Moxi GO II instrument (Orflo, Cat # MXG102). Cells were diluted to a final concentration of 3x10<sup>5</sup> cells/mL.

## Sample Preparation

For sample preparation, 1 mL of HepG2 cells was added to a 2 mL Hard Tissue Homogenizing Mix tube (OMNI, Cat # 19-628). The HepG2 samples were homogenized on the Bead Ruptor Elite™ Bead Mill Homogenizer (OMNI, Cat # 19-042E) at 4.5 m/s for 30 seconds. After homogenization, the lysate was centrifuged at 10,000 xg for 5 minutes to pellet cell debris.

## Workflow Integration

The Bead Ruptor Elite™ Bead Mill Homogenizer can process 24 or 48, 2 mL tubes at one time using either the 24 or 48-position 2 mL tube carriage (OMNI, Cat # 19-373) (OMNI, Cat # 19-378), respectively. After homogenization, the lysate is manually transferred into chemagen 24 or 96 well plates, depending on throughput. Alternatively, both 2 mL tube carriages have been formatted in-house with the capability to attach to a custom-designed support on the JANUS® G3 deck, along with Winprep integration as a custom Labware file. With this custom integration, automated liquid handling of lysate direct from the 24 or 48-position carriages into desired lysate plates is facilitated as a part of chemagic™ 360 automated nucleic acid extraction workflows.



**Figure 1.** Example Winprep deck layout showcasing the 24 position 2 mL tube carriage (top right, black) and 48 position 2 mL tube carriage (bottom left, pink).

## Nucleic Acid Extraction and Analysis

After centrifugation step, 100 µl of cell lysate was transferred to 12 separate wells of a 96-well plate (PerkinElmer, Cat # CMG-555-15) using the JANUS® G3 automated workstation and customized Reagent Liquid Transfer Winprep program. Similarly, 100 µl Lysis Buffer provided with the chemagic Tissue DNA Kit (PerkinElmer, Cat # CMG-723) was transferred to each well containing the lysate. Next, 6 µl of Proteinase-K, included in the chemagic Tissue DNA kit, was manually transferred to all wells along with 5 µl RNase A (Thermo Scientific, Cat # EN0531). The plate was then incubated at 37°C for 10 minutes. During this time, Elution Buffer and Magnetic Beads were transferred into respective plates using a Reagent Liquid Transfer Protocol on the JANUS® G3 automated workstation. After incubation, the tissue DNA extraction was carried out using the chemagic 360 nucleic acid extractor.


After completion of the chemagic™ 360 extraction protocol, eluted nucleic acid concentration and integrity was determined by  $A_{260}/A_{280}$  spectrophotometry.

Quantitative PCR

Ten microliters of eluted HepG2 DNA was added to qPCR mix (BioRad, Cat # 1725122) along with 5  $\mu$ M of forward and reverse 18S primers. The 18S gene was targeted with forward primer 5' - CAG CAG CCG CGG TAA TTC C – 3', reverse primer 5' - CCC GTG TTG AGT CAA ATT AAG C – 3' yielding a product size of 676 bp. An 18S-positive DNA extract was used as the positive control along with nuclease-free water for the negative control. Reactions were loaded into the BioRad CFX Connect Real Time Instrument (BioRad Cat. # 1855201), amplified for 40 cycles and visualized via gel electrophoresis.

Results

The outlined experiments using the chemagic™ Tissue DNA kit and Bead Ruptor Elite™ Bead Mill Homogenizer for sample preparation resulted in nucleic acid that is suitable for downstream PCR analysis. Eluted nucleic acid had an average yield of 3.9  $\mu$ g. Amplified HepG2 DNA had an average Cq value of 16.16, revealing that eluted DNA is intact and compatible with downstream analysis (Table 1).



**Figure 2.** HepG2 DNA concentration from 12 samples. Lysate from samples 1-12 was prepared with Pro-K and RNase A. Average  $A_{260}/A_{280}$  from samples 1-12 was 2.21.

Sample	Cq
HepG2 1	16.05
HepG2 2	16.32
HepG2 3	15.08
HepG2 4	16.04
HepG2 5	15.28
HepG2 6	15.75
HepG2 7	16.35
HepG2 8	16.49
HepG2 9	16.11
HepG2 10	16.01
HepG2 11	16.02
HepG2 12	18.46
Positive Extract Control	22.01
Negative Control	39.07

**Table 1.** 18S qPCR Cq values from extracted HepG2 DNA.

## Conclusions

In the enclosed application note, we have showcased proof-of-concept integration capabilities of the Bead Ruptor Elite™ Bead Mill Homogenizer and chemagic™ 360 automated nucleic acid extractor for isolation of DNA from *in-vitro* grown HepG2 cells. The Bead Ruptor Elite™ Bead Mill Homogenizer and 2 mL Hard Tissue Homogenizing Mix was demonstrated as a solution for sample preparation of *in-vitro* cultured cells, producing a homogenate suitable for downstream automated nucleic acid extraction. After a 70-minute extraction, eluted DNA was high-yield and intact, proving suitable for downstream analysis via qPCR. Additionally, the Bead Ruptor Elite™ Bead Mill Homogenizer can accommodate a variety of sample preparation throughput demands with either 24 or 48-position 2 mL tube carriages. Using either manual methods or custom-integrated JANUS® G3 liquid handling, homogenate transfer and extraction preparation is achieved.

## Ordering Information

Equipment	Catalog Number
Bead Ruptor Elite™ Bead Mill Homogenizer	19-042E
Bead Ruptor Elite™ Bead Mill Homogenizer 2 mL Tube Carriage	19-373
Bead Ruptor Elite™ Bead Mill Homogenizer Finger Plate	19-370
Bead Ruptor Elite™ Bead Mill Homogenizer 48 position 2 mL Tube Carriage	19-378
Bead Ruptor Elite™ Bead Mill Homogenizer 48 position 2 mL Tube Carriage Finger Plate	19-370-248
Hard Tissue Homogenizing Mix 2.8 mm Ceramic Beads	19-628

For research use only. Not for use in diagnostic procedures.

Contact Info for your local  
OMNI International, a PerkinElmer company  
Representative:  
[sales@OMNI-inc.com](mailto:sales@OMNI-inc.com) or  
<https://us.OMNI-inc.com/contact-us>

**PerkinElmer, Inc.**  
940 Winter Street  
Waltham, MA 02451 USA  
P: (800) 762-4000 or  
(+1) 203-925-4602  
[www.perkinelmer.com](http://www.perkinelmer.com)

**OMNI International, a PerkinElmer company**  
935 Cobb Place Blvd NW  
Kennesaw, GA 30144 USA  
P: (800) 766-4431 or  
(+1) 770-421-0058  
[www.omni-inc.com](http://www.omni-inc.com)



For a complete listing of our global offices, visit [www.perkinelmer.com/ContactUS](http://www.perkinelmer.com/ContactUS)

Copyright ©2021, PerkinElmer, Inc. All rights reserved. PerkinElmer® is a registered trademark of PerkinElmer, Inc. All other trademarks are the property of their respective owners.

12/20/2021