

# Omni Mixer Series

## HOMOGENIZERS



Homogenizing tough  
samples for over 60 years





## Mixer Homogenizer

ORDER NO. 17105 (115V)

ORDER NO. 17106 (220V)

**Processing Range:** 0.25 mL to 10 L

**Generator Probes:** 10, 20, 30 mm

**Blade Assemblies:** 1" and 2" in sealed chambers

**Power Rating:** 600 watt

**Speed:** Variable from 500 to 18,000 rpm

**Dimensions (with stand):** 24" H x 10" W

**Weight:** 9.5 kg (21 lbs.)

**Sound Level:** <68 db

**Warranty:** 3 year motor warranty



## Macro ES Homogenizer

ORDER NO. MES-115 (115V)

ORDER NO. MES-220 (220V)

**Processing Range:** 0.03 mL to 30 L

**Generator Probes:** 10, 20, 30, 45, 55 mm

**Blade Assemblies:** 1" and 2" in sealed chambers

**Power Rating:** 1,800 watt

**Speed:** Variable from 1,000 to 20,000 rpm, soft start, 0.5% accuracy

**Programs:** 12 programs with 12 time and speed steps, including ramping function

**Interface:** RS-232 and Windows software

**Dimensions:** 34.5" max

**Weight:** 14.5 kg (32 lbs.)

**Sound Level:** <68 db

**Warranty:** 3 year motor warranty

**Compliance:** CE Approved

## Macro Homogenizer

ORDER NO. 17505 (115V)

ORDER NO. 17506 (220V)

**Processing Range:** 0.25 mL to 30 L

**Generator Probes:** 10, 20, 30, 45, 55 mm

**Blade Assemblies:** 1" and 2" in sealed chambers

**Power Rating:** 1,800 watt

**Speed:** Variable from 1,000 to 20,000 rpm

**Dimensions (with stand):** 24" H x 10" W

**Weight:** 9.5 kg (21 lbs.)

**Sound Level:** <68 db

**Warranty:** 3 year motor warranty



935 Cobb Place Blvd.  
Kennesaw, GA 30144  
800.776.4431 • 770.421.0058  
www.omni-inc.com

REV.E