

Tissue Master

TISSUE HOMOGENIZER

The **Tissue Master 125** is an extremely economical homogenizer that is well suited for customers that are occasionally homogenizing small volume samples with little variation in sample size. The Tissue Master 125 comes complete with a 7 mm stainless steel probe.

The **Omni Tissue Master 125 Homogenizer Complete Kit** includes a 5 mm, 7 mm and 10 mm stainless steel generator probes, specifically designed for processing volumes of <math><200 \mu\text{L}</math> to 100 mL. These probes can handle a wide variety of different sample sizes, making the **Tissue Master** a great choice for a laboratory with many applications, but a limited budget.

Processing Range: <math><200 \mu\text{L}</math> to 100 mL

Speed: Variable from 500 to 35,000 rpm

Power Rating: 125 watt

Weight: 535 g (18.8 oz.)

Sound Level: <math><72 \text{ db}</math>

Warranty: 1 year motor warranty

Standards Approval: CE Approved

Features

- Powerful 125 watt motor
- Light weight ergonomic design
- Value priced, complete with 7 mm probe
- Complete kit includes 5 mm, 7 mm and 10 mm stainless steel generator probes



Omni Stainless Steel generator probes

OMNI TISSUE MASTER 125 HOMOGENIZER COMPLETE KIT

(WITH 5, 7 AND 10 MM PROBES) ALSO AVAILABLE INDIVIDUALLY

Available in 115V or 220V

ORDER NO. TM125-115

TM125-220



5 MM PROBE



7 MM PROBE



10 MM PROBE



VALUE PRICED HOMOGENIZING KIT

WITH 5 MM PROBE ORDER NO. TM125-115-05

WITH 7 MM PROBE ORDER NO. TM125-115

WITH 10 MM PROBE ORDER NO. TM125-115-10



OMNI
International

935-C Cobb Place Blvd
Kennesaw, GA 30144
Phone: 770.421.0058 • Fax: 770.421.0206
www.omni-inc.com



REV.B